

Radiolucent NeoLead®

ECG Electrodes

INDICATIONS FOR USE

The Radiolucent NeoLead and Micro NeoLead ECG Electrodes are intended for ECG monitoring. They are intended for use on pediatric (neonates, infants and children) patients.

DIRECTIONS FOR USE

Preparation

Clean and dry skin per hospital protocol.

Note: It is important for skin to be as clean and dry as possible. Electrodes will not adhere properly to moist skin or hair. Oils and lotions will also adversely affect adhesion.

DO NOT USE ALCOHOL.

Application

Step 1

Before removing liner, warm the lead in hands for 60 seconds. May also be held under radiant warmer or heel warmer for 10 to 15 seconds.

Step 2

Peel and discard liner.

Step 3

Place the lead on skin and hold in place for 60 seconds to ensure proper adhesion.

Removal

Step 1

Replace NeoLead every 5 days or per hospital protocol, whichever is sooner.

Step 2

Saturate the lead with water or saline.

Step 3

Slowly peel the lead away from the skin as you swab with water or saline.

GENERAL INFORMATION

Product Materials

Electrode: Silver silver chloride, hydrogel, NeoBond Hydrocolloid, non-woven fabric
Wire: Copper, polypropylene
Made with natural rubber latex: No
Made with plasticizer DEHP: No

Size

N305, N315
Diameter: 26 mm
Wire length: 61 cm

N306

Electrode width: 25 mm
Electrode height: 10 mm
Wire length: 61 cm

Packaging

Sets individually packaged: Yes
Package type: Pouch
Dispenser box: Yes

Visit neotechproducts.com for the most up-to-date information.

Pouches are made from paper, polyethylene and foil. Boxes are printed with environmentally friendly soy ink on recyclable material.

CAUTIONS & WARNINGS

Federal law restricts this device to sale by or on the order of a physician.

Radiolucent NeoLead and Micro NeoLead ECG Electrodes should only be used while patient is under continuous direct supervision of healthcare professionals.

Electrodes will not adhere properly to moist skin or hair. Oils and lotions will adversely affect adhesion.

Do not use alcohol to clean skin prior to application.

Do not use Radiolucent NeoLead (N305, N306, N315) during MR procedures.

Discontinue immediately if skin irritation occurs.

Disposable. Do not reuse.

Single patient use only.



For more information or to order your FREE sample call
1-800-966-0500 or visit neotechproducts.com

Neotech Products LLC • 28430 Witherspoon Parkway, Valencia, CA 91355