DIRECTIONS FOR USE Non-Sterile/Single-Patient-Use Only DO NOT CLEAN OR STERILIZE

INDICATIONS FOR USE: The intended use of the Neo-StatCO₂<Kg[®] is to provide a semi-quantitative visualization of the CO₂ in the patient airway.

DESCRIPTION: The Mercury Medical Neo-StatCO₂<Kg End-Tidal CO₂ Detector is a fast, durable colorimetric breath indicator for visualization of exhaled CO₂.

WARNINGS: 1. Do not use for the detection of hypercapnia. 2. Do not use for the detection of main-stem bronchial intubation. 3. Do not use during mouth-to-tube ventilation.


Diagram #1

PRECAUTIONS: 1. Single Use Only: Do not sterilize or reuse as the intended function will be affected. 2. Do not use unit if package is not sealed or after more than 24 hours exposure to ambient air.

- 1. Check to ensure that the package is intact. Check the expiration date located on the outer carton label and foil pack before use. 2. Insert the Endotracheal Tube (ET tube) into the patient and inflate the cuff.

DIRECTIONS FOR USE - PLEASE FOLLOW THESE DIRECTIONS CAREFULLY: 1. Check to ensure that the package is intact. Check the expiration date located on the outer carton label and foil pack before use.


- 1. Indicator Color: Blue No Carbon Dioxide is present 2. Green Exposed to 1.0 - 2.0% CO2 3. Yellow Green At intermediate CO2 concentrations of 2.0% - 5.0%.

WARNING: If, after six breaths, the Mercury Medical Neo-StatCO₂<Kg End-Tidal CO₂ Detector does not display cyclical color changes in phase with respiration, this is a strong indication of incorrect intubation.

Table with 2 columns: Indicator Color, Indicators. Lists color codes and corresponding CO2 concentrations.

If used for longer than twenty-four hours, the indicator colors eventually fade to a yellow-light green.

Table with 2 columns: MECHANICAL SPECIFICATIONS, Connector Ports. Details internal volume, resistance, and weight.

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Mini StatCO₂[®]

End Tidal CO₂ Detector for Patients
with body weight between 1 and 15 Kg.
(2.2 - 33 lb.)

(en)

DIRECTIONS FOR USE
Non-Sterile/Single-Patient-Use Only
DO NOT CLEAN OR STERILIZE

DESCRIPTION:
The Mercury Medical Mini StatCO₂[®] End-Tidal CO₂ Detector is a fast, durable colorimetric breath indicator for visualization of exhaled CO₂. The detector is designed to connect between an endotracheal tube and a breathing device to verify proper intubation. Exhaled gas passes through the indicator to detect approximate ranges of end-tidal CO₂ by color comparison. Color change is very rapid and the detector has the ability to show color changes up to a respiratory rate of 50 breaths per minute. The detector may be used during patient transport or in locations where intubations are performed. While the indicator will function at relative humidities (RH) of up to 100%, it is not recommended for use below 10% relative humidity. The Mercury Medical Mini StatCO₂[®] End-Tidal CO₂ Detector may be used for up to twenty-four hours.

INDICATIONS:
The intended use of the Mini StatCO₂ is to provide a semi-quantitative visualization of the CO₂ in the patient airway. It is an adjunct in patient assessment to be used in conjunction with other methods to determine clinical signs and symptoms by and on the order of a physician.

CONTRAINDICATIONS:
- Do not use for the detection of hypercapnia,
- Do not use for the detection of main-stem bronchial intubation,
- Do not use during mouth-to-tube ventilation,
- Do not use the CO₂ detector to detect oropharyngeal tube placement,
- Standard clinical assessment must be used.

When low pulmonary perfusion coincides with accidental esophageal intubation, colorimetric CO₂ indication cannot be properly interpreted. However, if proper tube placement is ascertained by independent means, then the detector may be used to help assess the progress of positive pressure ventilation as evidenced by an increase in end tidal CO₂.

WARNINGS: (1) The Mini StatCO₂ Detector cannot be interpreted properly by individuals who experience blue-yellow color blindness. (2) Do not use this device on patients with body weight greater than 15 kg (33 lb) due to possibility of increased airway resistance. (3) Use this device with air and oxygen only; chemical interactions will affect device accuracy. (4) Do not use this device with a nebulizer. (5) Before using the Mercury Medical Mini StatCO₂ End Tidal CO₂ Detector, please read entire contents of this Directions for Use. (6) Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DIRECTIONS FOR USE - PLEASE FOLLOW THESE DIRECTIONS CAREFULLY:

- Check to ensure that the package is intact. Check the expiration date located on the outer carton label and foil pack before use. Remove the device from the package and check the color of the detector by comparing it to the reference scale. Match the initial color of the indicator to the blue color label on the detector window. (See Diagram #1) If the blue color of the detector is not the same color as the blue label, or if it is lighter than the blue label, **DISCARD THE DEVICE AND DO NOT USE.**
- Insert the Endotracheal Tube (ET tube) into the patient and inflate the cuff (if equipped) per manufacturer's guidelines.
- Firmly attach the device between the ET tube and the breathing device. Pull the red tab from the device to activate before use. (See Diagram #2)
- Wait until six breaths of normal tidal volume have occurred. Initial Mini CO₂ Detector color (yellow) may be interpreted as a false positive if read before delivery of six breaths.
- After six breaths, compare the indicator color in the device window on full expiration to readings printed on the detector cover. During use, the detector will continue to change color in phase with respiration.

Diagram #1

Diagram #2

WARNING:
If, after six breaths, the Mercury Medical Mini StatCO₂[®] End-Tidal CO₂ Detector does not display cyclical color changes in phase with respiration, this is a strong indication of incorrect intubation. If it ceases to change color, this indicates an equipment malfunction (such as a disconnection, CO₂ rebreathing or a patient-related problem). Use appropriate alternative procedure.
6. Note that the detector color will continue to fluctuate from inspiration to expiration for up to twenty-four hours. When this gas contains CO₂ (i.e., expired gas from the patient), the indicator color will change. Following are the color indications of the Mini StatCO₂ indicator device.

Indicator Color:	Indicates
Blue	No Carbon Dioxide is present
Green	Exposed to 1.0 - 2.0% CO ₂
Yellow Green	At intermediate CO ₂ concentrations of 2.0% - 5.0%.
Yellow	Exposed to 5.0% or more CO ₂
Permanent Yellow	A damaged indicator

7. Remove and discard the Mercury Medical Mini StatCO₂ End-Tidal CO₂ Detector after use, not to exceed 24 hours.

WARNING:
Periodic color changes reflect the breathing pattern of the patient. However, a permanent blue to blue-green color indicates a lack of exhaled carbon dioxide, the cause of which requires immediate attention. If the blue-green color does not appear during a breathing cycle, this may indicate a significant degree of carbon dioxide rebreathing and should also be of immediate concern.

Limitations of CO₂ Detector Function:

- The detector responds almost instantaneously when exposed to CO₂. It can be used up to a respiratory rate of 50 breaths per minute.
- Concentrations as low as 1.0% - 2.0% will cause a noticeable change in color from blue to green. At concentrations of approximately 5%, the color will be yellow, and will remain yellow even if the CO₂ concentration should rise to higher levels. If used for longer than twenty-four hours, the detector color will eventually fade to a yellow-light green.
- Device should not be used by persons with blue/yellow color blindness.

CAUTION:

- Single Use Only: Do not sterilize or reuse as the intended function will be affected. This may increase the risk of possible transmission of infection. Discard after use.
- Do not use unit if package is not sealed or after more than 24 hours exposure to ambient air. If the Mercury Medical Mini StatCO₂ End-Tidal CO₂ Detector does not produce conclusive results after intubation, the ET tube should immediately be inserted unless correct anatomical location can be confirmed by other means.
- Do not use in the presence of acidic liquid or medication.
- Avoid exposure to strong sunlight or other sources of Ultraviolet light.
- Revent liquid water from covering indication surface, since this may impair the device's function.
- Store in a sealed bag at or below room temperature. Storing above room temperature for long durations may decrease shelf life.
- This unit does not take the place of traditional ET CO₂ monitoring which remains necessary to provide quantitative measurements and patient alarms.

MECHANICAL SPECIFICATIONS:	
Internal Volume: 3 cc Resistance to Flow: 2.5 cm H ₂ O nominal at 10 L/min flow Weight: 5 grams nominal	Connector Ports: Patient end: 15 mm I.D. standard taper Circuit end: 15 mm O.D. standard taper

CE MD Do Not Reuse Not Made with Natural Rubber Rx By Prescription Only Non Sterile Keep Dry Phthalates EC/REP EMERGO EUROPE Prinsesgracht 20 2514 AP The Hague The Netherlands Caution

MercuryMedical
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercury-med.com

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MercuryMedical **StatCO₂**[®] End Tidal CO₂ Detector for Patients with body weight more than 15 Kg. **(en)**

11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

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DESCRIPTION:
The Mercury Medical StatCO₂[®] End-Tidal CO₂ Detector is a fast, durable colorimetric breath indicator for visualization of exhaled CO₂. The detector is designed to connect between an endotracheal tube and a breathing device to help verify proper intubation. Exhaled gas passes through the indicator to detect approximate ranges of end-tidal CO₂ by color comparison. Color change is very rapid and the detector has the ability to show color changes up to a respiratory rate of 50 breaths per minute. The detector may be used during patient transport or in locations where intubations are performed. While the indicator will function at relative humidities (RH) of up to 100%, it is not recommended for use below 10% relative humidity. The Mercury Medical StatCO₂ End-Tidal CO₂ Detector may be used for up to twenty-four hours.

INDICATIONS:
The intended use of the StatCO₂ is to provide a semi-quantitative visualization of the CO₂ in the patient airway. It is an adjunct in patient assessment to be used in conjunction with other methods to determine clinical signs and symptoms by and on the order of a physician.

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
WARNINGS: (1) The StatCO₂ Detector cannot be interpreted properly by individuals who experience blue-yellow color blindness. (2) Do not use this device on patients with body weight less than 15 kg (33 lb) due to the risk of rebreathing exhaled CO₂. (3) Use this device with air and oxygen only; chemical interactions will affect device accuracy. (4) Do not use this device with a nebulizer. (5) Before using the Mercury Medical StatCO₂ End Tidal CO₂ Detector, please read the entire contents of this Directions for Use. (6) Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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- Check to ensure that the package is intact. Check the expiration date located on the outer carton label and foil pack before use. Remove the device from the package and check the color of the detector by comparing it to the reference scale. Match the initial color of the indicator to the blue color label on the detector window. (See Diagram #1) If the blue color of the detector is not the same color as the blue label, or if it is lighter than the blue label, **DISCARD THE DEVICE AND DO NOT USE.**
- Insert the Endotracheal Tube (ET tube) into the patient and inflate the cuff (if equipped) per manufacturer's guidelines.
- Firmly attach the device between the ET tube and the breathing device. Press the red tab from the device to activate before use. (See Diagram #2)
- Wait until six breaths of normal tidal volume have occurred. Initial CO₂ Detector color (yellow) may be interpreted as a false positive if read before delivery of six breaths.
- After six breaths, compare the indicator color in the detector window on full expiration to ratings printed on the detector cover. During use, the detector will continue to change color in phase with respiration.

WARNING:
If, after six breaths, the Mercury Medical StatCO₂[®] End-Tidal CO₂ Detector does not display cyclical color changes in phase with respiration, this is a strong indication of incorrect intubation. If it ceases to change color, this indicates an equipment malfunction (such as a disconnection, CO₂ rebreathing or a patient-related problem). Use appropriate alternative procedure.

6. Note that the detector color will continue to fluctuate from inspiration to expiration for up to twenty-four hours. When this gas contains CO₂ (i.e., expired gas from the patient), the indicator color will change. Following are the color indications of the StatCO₂ indicator device.

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7. Remove and discard the Mercury Medical StatCO₂ End-Tidal CO₂ Detector after use, not to exceed 24 hours.

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- The detector responds almost instantaneously when exposed to CO₂. It can be used up to a respiratory rate of 50 breaths per minute.
- Concentrations as low as 1.0% - 2.0% will cause a noticeable change in color from blue to green. At concentrations of approximately 5%, the color will be yellow, and will remain yellow even if the CO₂ concentration should rise to higher levels. If used for longer than twenty-four hours, the detector colors will eventually fade to a yellow-light green.
- Device should not be used by persons with blue/yellow color blindness.

CAUTION:
- Single Use Only: Do not sterilize or reuse as the intended function will be affected. This may increase the risk of possible transmission of infection. Discard after use.
- Do not use unit if package is not sealed or after more than 24 hours exposure to ambient air. If the Mercury Medical StatCO₂ End-Tidal CO₂ Detector does not produce conclusive results after intubation, the ET tube should immediately be reinserted unless correct anatomical location can be confirmed by other means.
- Do not use in the presence of acidic liquid or medication.
- Avoid exposure to strong sunlight or other sources of ultraviolet light.
- Prevent liquid water from covering indication surface, since this may impair the device's function.
- Store in a sealed bag at or below room temperature. Storing above room temperature for long durations may decrease shelf life.
- This unit does not take the place of traditional ET CO₂ monitoring which remains necessary to provide quantitative measurements and patient alarms.

MECHANICAL SPECIFICATIONS:	
Internal Volume: 25 cc	Connector Ports: Patient end: 15 mm I.D. standard taper
Resistance to Flow: 3.0 cm H ₂ O nominal at 60 L/min flow	Circuit end: 15 mm O.D. standard taper
Weight: 22 grams nominal	

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©Marchio di fabbrica USA della Mercury Enterprises, Inc. Brevetto n. US 7.578.971 B2, US 8.256.414 B2, US 8.402.964 B2.