



FLEXI-Board

Comfortable & Flexible
Waterproof Armboard

Flexi board is constructed from soft polyurethane foam with an aluminium strip core, they can be easily shaped and positioned on the patient without losing their profile.

FLEXI-board is available in three sizes: Small, Medium and Large.

Features

- Malleable soft waterproof foam cover.
- 3 different sizes available.
- Simple application and removal.
- Maximising comfort to limbs.
- Aluminium core ensures the profile is maintained.

Ordering Information

Cat No.	Sizes	Qty./Box
743510	Small	100
743511	Medium	80
743512	Large	60

RESPICARE SOLUTIONS