


LATEX FREE



Neo-FLOW[®] Oxygen Nasal Cannulas

Hand dipped Smoothflow cannulas
For Total Patient Comfort

Neo-FLOW Oxygen Nasal Cannulas

The NeoFLOW Oxygen Nasal cannula is an extremely soft cannula with Curved tips specially designed for total infant / child comfort.


Hand-dipped flexible contoured smooth flow design ensures the patient feels no discomfort whatsoever.

Kink resistant tubing with inner star lumen for uninterrupted flow of oxygen.

Over the ears lit secures position on patient.

Cannula Size	Tube Length	Product Ref.
Paediatric	2.1 metre	NCP 100
Infant	2.1 metre	NCI 100
Neonate	2.1 metre	NCN 100

Oxygen Cannula



Product Ref.	1	2			3			4			5
	Cannula Size	Nasal Prong "D"		Nasal Prong "L" Length	Cannula Tube# Length		Cannula Tube# I.D. O.D.		Oxygen Tubing* Length I.D. O.D.		Oxygen Connector Style
NCP 100	Paediatric	3.25	2.5	11.0 ± 1.0	450 ± 10	3.2 ± 0.1	1.9 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange
NCI 100	Infant	3.0	2.0	10.0 ± 1.0	450 ± 10	3.2 ± 0.1	1.9 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange
NCN 100	Neonate	2.0	1.5	5.5 ± 0.5	250 ± 10	2.2 ± 0.1	1.1 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange

All Dimensions in millimeters

- Soft PVC Tube * - Star Lumen Tube

RESPICARE SOLUTIONS

260, Adi Mansion, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai - 400 002.
Tel.No.: (022) 22061706, 22017197, 22005405, +91 98200 62869. Email : cncindia@vsnl.com • www.cncmedical.net

