


LATEX FREE



Neo-FLOW[®] Oxygen Nasal Cannulas

Hand dipped Smoothflow cannulas
For Total Patient Comfort

Neo-FLOW Oxygen Nasal Cannulas

The NeoFLOW Oxygen Nasal cannula is an extremely soft cannula with Curved tips specially designed for total infant / child comfort.


Hand-dipped flexible contoured smooth flow design ensures the patient feels no discomfort whatsoever.

Kink resistant tubing with inner star lumen for uninterrupted flow of oxygen.

Over the ears lit secures position on patient.

Tube Length	Cannula Size	Product Ref.
2.1 metre	Paediatric	NCP 100
2.1 metre	Infant	NCI 100
2.1 metre	Neonate	NCN 100

Oxygen Cannula



Product Ref.	1 Cannula Size	2 Nasal Prong			3 Cannula Tube [#]			4 Oxygen Tubing [*]			5 Oxygen Connector Style
		"D" I.D.	"D" O.D.	"L" Length	Length	I.D.	O.D.	Length	I.D.	O.D.	
NCP 100	Paediatric	3.25	2.5	11.0 ± 1.0	450 ± 10	3.2 ± 0.1	1.9 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange
NCI 100	Infant	3.0	2.0	10.0 ± 1.0	450 ± 10	3.2 ± 0.1	1.9 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange
NCN 100	Neonate	2.0	1.5	5.5 ± 0.5	250 ± 10	2.2 ± 0.1	1.1 ± 0.1	2100 ± 50	6.3 ± 0.1	4.3 ± 0.1	5-8mm Flange

All Dimensions in millimeters

- Soft PVC Tube * - Star Lumen Tube

RESPICARE SOLUTIONS

260, Adi Mansion, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai - 400 002.
Tel.No.: (022) 22061706, 22017197, 22005405, +91 98200 62869. Email : cncindia@vsnl.com • www.cncmedical.net

