



RESPI – EVAC[®]

Endotracheal Tube

With integrated suction lumen



Product Introduction

Ventilator Associated Pneumonia (VAP) is one of the most common nosocomial infection that develops more than 48 hours after endotracheal intubation. Clinical Studies have shown that one of the most common cause is contamination of the lower respiratory tract due to seepage of the secretions accumulated above the Endotracheal cuff.

The RESPI-EVAC[®] VAP Control Strategy

Continuous Aspiration of Subglottic Secretions (CASS) is increasingly gaining acceptance as an effective clinical strategy to prevent the onset of VAP.

The Respi-Evac[®] Endotracheal tube with integrated suction lumen and evacuation port is a simple, convenient and cost effective solution which enables clinicians to manage secretions above the cuff by

- Removing secretions that may collect above the cuff.
- Reducing the risk of secretions moving past the cuff.
- Preparing the patient for extubation.

Product Benefits

- Large evacuation hole : Close to the cuff to help removal of secretions.
- Yellow suction connector: Easy to identify/good grip to attach to suction tubing.
- Yellow tethered dust-cap : Prevents contaminants from entering when not in use.
- Clear suction line : Allows visualisation of secretions.
- Cuff : High volume low pressure cuff
: Thin cuff creates low pressure seal.

At a Glance

- Endotracheal tube with integrated suction lumen.
- Effective Strategy for managing VAP by enabling CASS.
- Integral suction lumen allows continuous suctioning without risking trauma to the vocal cords as with manual catheter suctioning.
- Enables reduction in patient hospital stay & reduction in patient care cost associated with VAP.
- Manufactured in a USFDA approved facility.



Item code	Size id
700600	6.0mm
700650	6.5mm
700700	7.0mm
700750	7.5mm
700800	8.0mm
700850	8.5mm
700900	9.0mm

RESPICARE SOLUTIONS

260, Adi Mansion, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai - 400 002.
Tel.No.: (022) 22061706, 22017197, 22005405, +91 98200 62869. Email : enquiry@respicare.co.in • www.respicare.co.in

