



RESPI - CUFF

Cuff pressure indicator
with graduation

Representing a significant advance in tracheal cuff pressure monitoring, the Respi-Cuff cuff pressure indicator is one of the most innovative intubation safety devices.

- ✓ Continuous visual display combined with graduation eliminates the need for constant spot-checking and reduces the risk of Post Ventilatory Pneumonia
- ✓ Negative pressure graduation ensures cuff is fully deflated, reducing the risk of trauma during removal
- ✓ Lightweight device reduces drag and potential tube displacement
- ✓ Suitable for use on a wide range of airway devices extending usage into many clinical areas
- ✓ Simple design for ease of use
- ✓ Pressure can be quickly and easily adjusted avoiding excessive cuff inflation

Ordering Information:



Product Code	Pack Size
700002	50

RESPI CARE SOLUTIONS

Pravin Niwas 271-279, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai 400002.

Tel No.: 022 22061706, 22005405, 22017197, +929820062869. Email: enquiry@respicare.co.in . www.respicare.co.in