



# RESPI – FLEX<sup>®</sup>

## Reinforced Endotracheal Tube



Designed to reduce the risk of kinking when patient head is in extended or flexed position

### RESPI – FLEX<sup>®</sup> Reinforced Endotracheal Tube

- Uniform Spiral Stainless Steel Reinforcing wire within tube wall to reduce risk of kinking.
- Full Length radio-opaque line to assess the tube location.
- 15mm connector bonded to the tube eliminates risk of accidental disconnection.
- Suitable for both oral & nasal intubation.
- Soft, thin walled cuff ensures effective and a traumatic seal.



Item code	Size id
701300	3.0 mm
701350	3.5 mm
701400	4.0 mm
701450	4.5 mm
701500	5.0 mm
701550	5.5 mm
701600	6.0 mm
701650	6.5 mm
701700	7.0 mm
701750	7.5 mm
701800	8.0 mm
701850	8.5 mm
701900	9.0 mm
701950	9.5 mm



### RESPI CARE SOLUTIONS

