



# Neo PAP<sup>®</sup>

## Infant Bubble CPAP System

Complete system for delivery of Bubble CPAP to spontaneously breathing patients.

### Features

- NeoPAP<sup>®</sup> Bubble CPAP system generates natural pressure oscillations by creating bubbles in a water column.
- Combined with optimum humidification, NeoPAP<sup>®</sup> provides gentle yet effective respiratory support to spontaneously breathing patients.
- Available as a Complete Kit including:
  - NeoPAP<sup>®</sup> Bubble CPAP Generator,
  - Safety Valve Kit,
  - NeoVENT<sup>®</sup> Heated wire Breathing Circuit,
  - Humidification Chamber.
  - NeoPAP-S Nasal CPAP Prongs.
- Cost effective, as kit can be used as both Bubble CPAP circuit as well as Ventilator circuit

### RESPI CARE SOLUTIONS

Pravin Niwas, 271/279, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai - 400 002.

Tel.No.: (022) 22061706, 22017197, 22005405, +91 98200 62869. Email : enquiry@respicare.co.in • www.respicare.co.in

# ORDER INFORMATION

**NeoPAP®  
CPAP Generator**




**NeoPAP®  
Safety Valve**



**NeoPAP-S Nasal CPAP Kit**



## NeoVENT® Disposable Neonatal Heated Wire Breathing Circuit



CAT No.	Product	Product Contents
701161	<b>NeoPAP® Bubble CPAP System</b>	<ul style="list-style-type: none"> <li>• NeoPAP® Bubble CPAP Generator</li> <li>• Safety Valve Kit</li> <li>• NeoVent® Single Heated wire Breathing Circuit</li> <li>• Humidification Chamber</li> </ul>
701225	<b>NeoVENT® Single</b>	Single Heated wire circuit with chamber Circuit configuration: 110x1, 30x1, 60x3
702100	<b>NeoPAP-S Nasal CPAP Kit Size 0</b>	Prongs for < 700 g, Bonnet
702101	<b>NeoPAP-S Nasal CPAP Kit Size 1</b>	Prongs for 700 – 1250g, Bonnet
702102	<b>NeoPAP-S Nasal CPAP Kit Size 2</b>	Prongs for 1250 – 2000g, Bonnet
702103	<b>NeoPAP-S Nasal CPAP Kit Size 3</b>	Prongs for > 2000g, Bonnet

## RESPICARE SOLUTIONS

Pravin Niwas, 271/279, Dr. Cawasji Hormusji Street, Marine Lines, Mumbai - 400 002.

Tel.No.: (022) 22061706, 22017197, 22005405, +91 98200 62869. Email : enquiry@respicare.co.in • www.respicare.co.in